

HIGH RESOLUTION FIELD EMISSION SCANNING ELECTRON MICROSCOPE: HITACHI S-4800

BEAM ALIGNMENT PROCEDURE

1. Set V_{acc} and I_e to desired value.
2. Move stage to desired working distance with **Z** manual knob on sample chamber. The distance is in mm; this is the distance between the pole piece and the sample surface.
3. Under the **SEM** tab in the software interface, set **WD** (working distance).
4. Select focus mode **UHR**.
5. Press **F2** on the keyboard while the cursor is over the view window to degauss the objective lens. This should be done every time focus is greatly changed (by changing working distance), or V_{acc} or I_e are changed.
6. Adjust **FOCUS/BRIGHTNESS/CONTRAST** knobs to obtain the best image possible.
7. Check that I_e has not dropped from selected value. If it has, press **SET**.
8. Click **Align** button along top row of screen to open align dialog box.

Note: In general you want to align the beam at twice the magnification that you will be using for your images.

9. Align beam:
 - a. Click the **Beam Align** radio button.
 - b. Adjust **BRIGHTNESS/CONTRAST** knobs to obtain a clear disc. Use **STIGMA/ALIGNMENT** knobs **X** and **Y** to center disc on the target.
10. Align aperture:
 - a. Click the **Aperture Align** radio button.
 - b. Use **STIGMA/ALIGNMENT** knobs **X** and **Y** to minimize motion in image.
11. Align Stigma Align.X and Stigma Align.Y:
 - a. Click the **Stigma Align.X** radio button.
 - b. Use **STIGMA/ALIGNMENT** knobs **X** and **Y** to minimize motion in image.
 - c. Repeat for **Stigma Align.Y** radio button.
12. Select **Off** radio button to turn off alignment functions.

13. Adjust **FOCUS** knobs for best image.

14. Correct lens astigmatism by adjusting **STIGMA/ALIGNMENT** knobs **X** and **Y** knobs for best image.

When the astigmatism is out of adjustment, the image will be “smeared”. Round objects will appear oblong. Adjust **STIGMA/ALIGNMENT** knobs **X** and **Y** so that this “smearing” is minimized.

15. Repeat (14) and (15) until best image is obtained.

16. In some cases the initial focusing done in (7) would not have produced a quality image suitable for the alignment in (9) to (13). If this was the case, after (16) go back to (9) and repeat the process.